



**POLYPEPTIDE (PP) SOURCES, SUPPLIES,
AND INJECTION INSTRUCTIONS**

Companies that provide PPS are new and more are being developed. Tennant Foundation has no financial interests in any maker or distributor of PPs, but we support and salute the companies involved with the new emergence of PPs.

A. POLYPEPTIDE SOURCES

1. Integrative Peptides, LLC
103 W. Main St.
Hudson, MI 49247
Ph: 443-732-5272
Email: cs@integrativepeptides.com
Website: integrativepeptides.com
2. Biotech Peptides
San Diego, CA
Ph: 619-453-0156
Email: support@biotechpeptides.com
3. AAPPTEC
Ph: 502-968-2223
Fax: 502-968-3338
Email: sales@aapptec.com
Website: peptide.com
4. Limitless Like Nootrophics
help@limitlesslifenuootrophics.com
5. Super Smart
Us.supersmart.com
6. Pharm Grade Peptides
Gradepeptides.com
7. Peptide Sciences
231 St. Rose Pkwy, Suite 62
Henderson, NV 89052
800-986-6401
8. Biotics Research
Clinicalnutritioncenters.com
724-591-8126
9. Paradigm Peptides
1511 N. Covent
Boulonnais, IL 60914
800-706-3449

B. SUPPLIES FOR SUBCUTANEOUS INJECTIONS

- a. Insulin or allergy syringes with 26 to 30 gauge needle
- b. Sterile or bacteriostatic water
- c. Alcohol wipes

Supplies for subcutaneous injections can be obtained from local pharmacies or over the internet.

C. PREPARATION FOR INJECTIONS

- a. Inject 3 cc (ml) of sterile or bacteriostatic water into the vial which contains dry powder. It is ok to totally fill the vial.
- b. Shake well to mix.
- c. Wipe the skin area to be injected with alcohol before injection.
- d. Initial injection is .10 to .20 cc (ml) subcutaneous with an insulin syringe.

D. STORAGE

Keep in a dark, cool location. May store in refrigerator. Do not freeze.